



Surgical Disinfectant

USES AND DILUTIONS GUIDE



SURGICAL ROOM



SURGICAL EQUIPMENT



SURGICAL BED

Pre Operative	Post Operative
	
1:1 up to 1:10 1:10 = 1 Klens to 1 cup of	1:20 Dilution 1:20 = 1 Klens to 4 cups

APPLICATIONS

1. To Clean all Instruments used for both Pre and Post operative sessions.
2. Helps to protect all SS Instruments from Corrosion and infection.
3. When used regularly can Disinfect the patient post operative Beds

DID YOU KNOW

KlenS can save you money! When diluted with water, KlenS is ready to tackle micro organisms in and around your application area. It is very much user friendly and will not harm hands, and any type of skin when used continuously.

Manufactured by :



Marketed by : SBR Power Advantage India Pvt Ltd
Chennai-116 , Tamil Nadu , India , Email Id : info@advantagemep.com

K-8, 8th Street, Ambattur Industrial Estate, Ambattur, Chennai-600 058.

NOTE: Dilutions are approximate. "Do not use on Plastic Instrument Panels, Leather, Suede, Unfinished Wood, Opals or Pearls. For best results, mix dilutions in a clean empty spray bottle or bucket. Spot test for colorfastness on carpet and fabrics.